


Ideal thermal properties

Temperature increase during activation




BOWA EASY

The low temperature solution
BOWA neutral electrodes offer maximum patient safety with low temperature increase during activation.

Thanks to a thin hydrogel layer, they ensure faster heat removal and a much lower temperature increase compared to the current market standard.



Conventional neutral electrodes




Literature
Edrich J., Cookson C. C.,
Electrosurgical dispersive electrodes
heat cutaneous and subcutaneous skin
layers, Medical Instrumentation, Vol 21,
p 81-86, 1987

McAdams E., Bioelectrodes, in
Encyclopedia of Medical Devices and
Instrumentation, John & Sons, Inc.,
p 120-166, 2006

Schröder D., Wetzel D.,
Vorkommismeldungen zur
Hochfrequenz-Chirurgie, MTD Medizin-
Technischer Dialog, Vol 12/2009,
p 54-56, 2009

Connecting cable



Plug, generator side	Plug, neutral electrode side	REF	For generator type
		380-050 (4.5 m)	BOWA Erbe ICC/VIO International Valleylab REM Conmed
		386-050 (4.5 m)	BOWA Erbe ICC/VIO International Valleylab NON-REM Conmed
		101-150 (4.5 m)	ALSA Eschmann
		294-050 (4.5 m)	Erbe VIO/ICC/ACC Erbe T-Serie Martin

BOWA
EINFACH SICHER

BOWA-electronic GmbH & Co. KG
Heinrich-Hertz-Strasse 4 – 10
72810 Gomaringen | Germany

Phone +49 (0) 7072-6002-0
Fax +49 (0) 7072-6002-33
info@bowa.de | bowa-medical.com



BOWA-BRO-11188-NEUTRAL-ELECTRODES-EN-2017-07 Printed in Germany
We reserve the right to make technical and structural amendments Copyright of BOWA-electronic, Gomaringen

NEUTRAL ELECTRODES
INNOVATIVE, SAFE, ECONOMICAL



BOWA
EINFACH SICHER


Type „EASY Standard“



Type „EASY Universal“



Type overview – one for every need



		SPLIT	NON-SPLIT	REM SPLIT	NON-REM NON-SPLIT
Universal	> 5 kg	REF without cable (100 pcs.)		REF with cable 3 m (15 pcs.)	
	EASY Universal A _{cond1} = 90 cm ² SPLIT	816-092		818-092	
	EASY Standard A _{cond1} = 110 cm ² SPLIT	816-112		818-112	
	A _{cond1} = 110 cm ² NON-SPLIT	815-110		817-110	
Adults		> 15 kg			
	EASY A _{cond1} = 160 cm ² SPLIT	816-162		818-162	
	A _{cond1} = 160 cm ² SPLIT	816-161		818-161	
	A _{cond1} = 140 cm ² SPLIT	816-141		818-141	
	A _{cond1} = 140 cm ² NON-SPLIT	815-140		817-140	
Children		5 kg - 15 kg			
	A _{cond1} = 70 cm ² SPLIT	816-071		818-071	
	A _{cond1} = 70 cm ² NON-SPLIT	815-070		817-070	
Newborns		< 5 kg			
	EASY A _{cond1} = 40 cm ² SPLIT	816-042		818-042	
	A _{cond1} = 40 cm ² NON-SPLIT	815-040		817-040	
Packing					

* Conductive surface
Compatible with all major electrosurgical units.