

# EXACTO™ - MUNDO

Cold Snare




Reduce the Risk of  
**INCOMPLETE  
RESECTION**

# EXACTO™ - MUNDO


## Cold Snare

The **EXACTO - MUNDO Cold Snare** is a 14mm wide dedicated cold snare offering a 89% larger snare area than the market leading cold snare.<sup>1</sup> This larger cold snare offers an option for piecemeal cold snare polypectomy or cold EMR of large polyps.

### KEY FEATURES & BENEFITS:




- EXACTO™ - MUNDO Cold Snare
- EXACTO™ Cold Snare



**44%** WIDER SNARE LOOP

**1.5X** LARGER SNARE SURFACE

## A larger snare area **REDUCES THE RISK OF INCOMPLETE RESECTION.**<sup>2</sup>



## Product Information Chart

PRODUCT NUMBER	DESCRIPTION	SHEATH DIAMETER (MM)	LENGTH (CM)	SNARE SIZE (MM)	UNITS/BOX
BX00711114	EXACTO™ - MUNDO Cold Snare	2.5	230	14	10
BX00711078	EXACTO™ - MUNDO Cold Snare	2.5	230	14	20
BX00711079	EXACTO™ - MUNDO Cold Snare	2.5	230	14	50

1. Data on file at STERIS.

2. Cronin et al. (2024). Geometry of cold snare polypectomy and risk of incomplete resection. *Endoscopy*, 56(3), 214-219. <https://doi.org/10.1055/a-2184-1609>