

INSUFFLATION TUBING KITS & PARTICULATE FILTERS



- **Simple and easy** way to keep your patient and valuable equipment safe
- Feel **confident** that the gas you are delivering to your patient's abdomen is **free from debris and contaminants**.
- **Eliminate cross contamination** between patients and prevent the backflow of blood and body fluids entering the insufflator

Product Code	Description	Box Quantity
IFT0001AD	Insufflation set w/ tubing, filter & male luer both ends	25
IFT002AD	Insufflation Filter with 15mm Wolf Connector	25
IFT0015AD	Insufflation Filter with Push fit connector	25

FEATURES

- Cost effective
- Sterile – Single use – Easy peel pouch
- Bi Directional Bacterial/Viral Filter – 99.9999% efficiency at 0.027 micron
- Excellent Performance – High flow rate – Low Resistance – Enables rapid inflation
- Flexible 3m tubing
- Luer Lock Connector, Wolf Connector or push fit connector
- Latex Free